


Product Name: Recombinant Human CCDC141 Protein
Catalog #: PEH7101



Summary

Name	CC141_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZP82
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QLELLEEDSM KWRAKAEEYG RTLRSVEYC AMRDEINELK DSFKDIKKKF NNLKFNYYTKK NEKSRNLKAL KYQIQQVDMY AEKMQALKRK MEKVSNTSD SFLNYPSDKV NVLLEVMKDL QKHVDDFDKV VTDYKKNLDL TEHFQEVIEE CHFWYEDASA TVVRVGKYST ECKTKEAVKI LHQQFNKFIA PSVPQQEERI QEATDLAQHL YGLEEGQKYI EKIVTKHKEV LESVTELCES LTELKELKQ GDVLKMNPNL EDFHYDYIDL LKEPAKNKQT IFNEERNKGQ VQVADLLGIN GTGEERLPQD LKVSTDKEGG VQDLLLPEDM LSGEYECVS PDDISLPPLP GSPESPLAPS DMEVEEPVSS SLSLHISSYG VQAGTSSPGD AQESVLPPP AFADACNDKR ETFSSHFERP YLQFKAEPPL TSGFVKEST

image



Background

Product Name: Recombinant Human CCDC141 Protein
Catalog #: PEH7101



Alternative Names

Coiled-coil domain-containing protein 141 (Coiled-coil protein associated with myosin II and DISC1)

Note

For Research Use Only , Not for Diagnostic Use