


Product Name: Recombinant Human ATP6V0D2 Protein
Catalog #: PEH7106



Summary

Name	VA0D2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N8Y2
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MLEGAELYFN VDHGYLEGLV RGCKASLLTQ QDYINLVQCE TLEDLKIHLQ TTDYGNFLAN HTNPLTVSKI DTEMKRRLCG EFEYFRNHSL EPLSTFLTYM TCSYMIDNVI LLMNGALQKK SVKEILGKCH PLGRFTEMEA VNIAETPSDL FNAILIETPL APFFQDCMSE NALDELNIEL LRNKLYKSYL EAFYKFKCNH GDVTAEVMCP ILEFEADRRA FIITLNSFGT ELSKEDRETL YPTFGKLYPE GLRLLAQAED FDQMKNVADH YGVYKPLFEA VGGSGGKTLE DVFYEREVQM NVLAFNRQFH YGVFYAYVKL KEQEIRNIVW IAECISQRHR TKINSYIPIL

image



Background

Alternative Names	V-type proton ATPase subunit d 2 (V-ATPase subunit d 2) (Vacuolar proton pump subunit d 2)
--------------------------	--



Product Name: Recombinant Human ATP6V0D2 Protein
Catalog #: PEH7106

Note

For Research Use Only , Not for Diagnostic Use