


Product Name: Recombinant Human SLC13A1 Protein
Catalog #: PEH7107



Summary

Name	S13A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BZW2
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AEAEEVATQM TYFNGSTNHG LEIDESVNGH EINERKEKTK PVPGYNNDTG KISSKVELEK NSGMRTKYRT

image



Background

Alternative Names Solute carrier family 13 member 1 (Renal sodium/sulfate cotransporter) (Na⁺)/sulfate cotransporter) (hNaSi-1)

Note

For Research Use Only , Not for Diagnostic Use