


Product Name: Recombinant Human TMEM72 Protein
Catalog #: PEH7108



Summary

Name	TMM72_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A0PK05
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EVLASPEQYT DPSSSAVSTT GSGDTEQTYT FHGALKEGPS SLFIHMKSL KGTKKPSALQ PPNTLMELSL EPADSLAKKK QVHFEDNLVR IVPSLAEGLD DGDSEPEETT SDTTPPIPPP

image



Background

Alternative Names Transmembrane protein 72 (Kidney-specific secretory protein of 37 kDa)

Note

For Research Use Only , Not for Diagnostic Use