


Product Name: Recombinant Human ERP27 Protein
Catalog #: PEH7111



Summary

Name	ERP27_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96DN0
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SDGPGAAQEP TWLTDVPAAM EFIAATEVAV IGFFQDLEIP AVPI LHSMVQ KFPGVSFSGIS TDSEVLTHYN ITGNTICLFR LVDNEQLNLE DEDIESIDAT KLSRFIEINS LH MVTEYNPV TVIGLFNSVI QIHLLIMNK ASPEYEENMH RYQKAAKLFQ GKILFILVDS GMKENGKVIS FFKLKESQLP ALAIYQTLDD EWDTLPTAEV SVEHVQN FCD GFLSGKLLKE

image



Background

Alternative Names	Endoplasmic reticulum resident protein 27 (ER protein 27) (ERp27) (Inactive protein disulfide-isomerase 27)
--------------------------	---

Note

Product Name: Recombinant Human ERP27 Protein
Catalog #: PEH7111



For Research Use Only , Not for Diagnostic Use