


Product Name: Recombinant Human MYH1 Protein
Catalog #: PEH7118



Summary

Name	MYH1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P12882
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKHATENKVK NLTEEMAGLD ETIAKLTKEK KALQEAHQQT LDDLQAEEDK VNTLTAKAIK LEQQVDDLEG SLEQEKKIRM DLERAKRKLE GDLKLAQEST MDIENDKQQL DEKLKKKEFE MSGLQSKIED EQALGMQLQK KIKELQARIE ELEEIEAER ASRAKAQQR SDLSRELEEI SERLEEAGGA TSAQIEMNKK REAEFQKMRR

image



Background

Alternative Names	Myosin-1 (Myosin heavy chain 1) (Myosin heavy chain 2x) (MyHC-2x) (Myosin heavy chain IIx/d) (MyHC-IIx/d) (Myosin heavy chain, skeletal muscle, adult 1)
--------------------------	--

Product Name: Recombinant Human MYH1 Protein
Catalog #: PEH7118



Note

For Research Use Only , Not for Diagnostic Use