



Product Name: Recombinant Human KLHL40 Protein
Catalog #: PEH7119



Summary

Name	KLH40_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q2TBA0
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FMISEEGAVA YDPAANECYC ASLSNQVPKN HVSLVTKENQ VFVAGGLFYN EDNKEDPMSA YFLQFDHLDSEWLGMPPLPS PRCLFGLGEA LNSIYVVGGR EIKDGERCLD SVMCYDRLSF KWGESDPLPY VVYGHTVLSH MDLVYVIGGK GSDRKCLNKM CVYDPKKFEW KELAPMQTAR SLFGATVHDG RIIVAAGVTD TGLTSSAEVY SITDNKWAPF EAFPQERSSL SLVSLVGTLY AIGGFATLET ESGELVPTEL NDIWRYNEEE KKWEGVLREI AYAAGATFLP VRLNVLCCLK

image



Background

Alternative Names	Kelch-like protein 40 (Kelch repeat and BTB domain-containing protein 5) (Sarcosynapsin)
--------------------------	--

Product Name: Recombinant Human KLHL40 Protein
Catalog #: PEH7119



Note

For Research Use Only , Not for Diagnostic Use