


Product Name: Recombinant Human RYR1 Protein
Catalog #: PEH7122



Summary

Name	RYR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P21817
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PQAGGEAQPA RAENEKDATT EKNKKRGFLF KAKKVAMMTQ PPATPTLPRL PHDVVPADNR DDPEIILNTT TYYYSRVVFA GQEPSCVWAG WVTPDYHQHD MSFDLSKVRV VVTVMGDEQG NVHSSLKCSN CYMVVWGGDFV SPGQQGRISH TDLVIGCLVD LATGLMTFTA NGKESNTFFQ VEPNTKLFPA VFVLPTHQNV IQFELGKQKN

image



Background

Alternative Names	Ryanodine receptor 1 (RYR-1) (RyR1) (Skeletal muscle calcium release channel) (Skeletal muscle ryanodine receptor) (Skeletal muscle-type ryanodine receptor) (Type 1 ryanodine receptor)
--------------------------	--

Product Name: Recombinant Human RYR1 Protein
Catalog #: PEH7122



Note

For Research Use Only , Not for Diagnostic Use