

Product Name: Recombinant Human MYOT Protein
Catalog #: PEH7124



Summary

Name	MYOTI_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBF9
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LPSQPDYNSS KIPSAMDSNY QQSSAGQPIN AKPSQTANAK PIPRTPDHEI QGSKEALIQD LERKCLKCKDT LLHNGNQRLT YEEKMARRLL GPQNAAAVFQ AQDDSGAQDS QQHNSEHARL QVPTSQVRSR

image



Background

Alternative Names Myotilin (57 kDa cytoskeletal protein) (Myofibrillar titin-like Ig domains protein) (Titin immunoglobulin domain protein)

Note

For Research Use Only , Not for Diagnostic Use