


Product Name: Recombinant Human MYH2 Protein
Catalog #: PEH7125



Summary

Name	MYH2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKX2
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KKKMEGDLNE MEIQLNHANR MAAEALRNYR NTQGILKDTQ IHLDDALRSQ EDLKEQLAMV ERRANLLQAE IEELRATLEQ TERSRKIAEQ ELLDASERVQ LLHTQNTSLI NTKKKLETDI SQMQGEMEDI LQEARNAEEK AKKAITDAAM MAEELKKEQD TSAHLERMKK NMEQTVKDLQ LRLDEAEQLA LKGGKKQIQK LEARVRELEG EVESEQKRNA EAVKGLRKHE RRVKELTYQT EEDRKNILRL QDLVDKLQAK VKSYKROAEE AEEQSNTNLA KFRKLQHELE EAEERADIAE SQVNKLRVKS REVHTKVISE

image



Background

Alternative Names Myosin-2 (Myosin heavy chain 2) (Myosin heavy chain 2a) (MyHC-2a) (Myosin heavy chain IIa) (MyHC-IIa) (Myosin heavy chain, skeletal muscle,

Product Name: Recombinant Human MYH2 Protein
Catalog #: PEH7125



adult 2)

Note

For Research Use Only , Not for Diagnostic Use