


Product Name: Recombinant Human DUSP29 Protein
Catalog #: PEH7133



Summary

Name	DUS29_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q68J44
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MTSGEVKTSL KNAVSSAKRL SPKMEEEGEE EDYCTPGAFE LERLFWKGSP QYTHVNEVWP KLYIGDEATA LDRYRLQKAG FTHVLNAAHG RWNVDTGPDY YRDMDIQYHG VEADDLPTFD LSVFFYPAAA FIDRALSDDH SKILVHCVMG RSRSATLVLA YLMIHKDRTL VDAIQQVAKN RCVLPNRFGL KQLRELDKQL VQQRSSQRQ DGEEEDGREL

image



Background

Alternative Names Dual specificity phosphatase 29 (Dual specificity phosphatase 27) (Dual specificity phosphatase DUPD1) (EC 3.1.3.16) (EC 3.1.3.48)

Note

Product Name: Recombinant Human DUSP29 Protein
Catalog #: PEH7133



For Research Use Only , Not for Diagnostic Use