



Product Name: Recombinant Human NLRP6 Protein
Catalog #: PEH7136



Summary

Name	NLRP6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P59044
Host	E.coli
Species	Human
Molecular Mass	64 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TGPR LAVARE LLLA ALEELS QEQLKRFRHK LRDVGPDGRS IPWGRLERAD AVDLAEQLAQ FYGPEPALEV ARKTLKRADA RDVAAQLQER RLQRLGLGSG TLLSVSEYKK KYREHVLQLH ARVKERNARS VKITKRFTKL LIAPESAAPE EAMGPAAEPE PGRARRSDTH TFNRLFRRDE

image



Background

Alternative Names	NACHT, LRR and PYD domains-containing protein 6 (Angiotensin II/vasopressin receptor) (PYRIN-containing APAF1-like protein 5)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NLRP6 Protein
Catalog #: PEH7136

