


Product Name: Recombinant Human NLRP6 Protein
Catalog #: PEH7137



Summary

Name	NLRP6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P59044
Host	E.coli
Species	Human
Molecular Mass	91 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QGPAGIGKTM AAKKILYDWA AGKLYQGQVD FAFFMPCGEL LERPGTRSLA DLILDQCPDR GAPVPQMLAQ PQRLIFLDG ADELPA LGGP EAAPCTDPFE AASGARVLGG LLSKALLPTA LLLVTTRAAA PGRLQGRLCS PQCAEVRGFS DKDKKKYFYK YFRDERRAER AYRFVKENET LFALCFVPFV CWIVCTVLRQ QLELGRDLSR TSKTTTSVYL LFITSVLSSA PVADGPRLQG DLRNLCRLAR EGVLGRRRAQF AEKELEQLEL RGSKVQTLFL SKKELPGVLE TEVTYQFIDQ SFQEFLAALS YLLEDGGVPR TAAGGVGTLL RGDAQPHSHL VLTTRFLFGL LSAERMRDIE RHFGCMVSER VKQEALRWVQ GQGQGC PGVA PEVTEGAKGL EDTEEP EEEEE EGEEPNYPLE

image



Background

Product Name: Recombinant Human NLRP6 Protein
Catalog #: PEH7137



Alternative Names

NACHT, LRR and PYD domains-containing protein 6 (Angiotensin II/vasopressin receptor) (PYRIN-containing APAF1-like protein 5)

Note

For Research Use Only , Not for Diagnostic Use