


Product Name: Recombinant Human AMER2 Protein
Catalog #: PEH7138



Summary

Name	AMER2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N7J2
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EAGPSCDKHV PGP GK PALS K KNP GVVAYQG GGEEMASPDE VDDTYLQEFW DMLSQTEEQG PEPQEGAAKV AAALETKVVP ETPKDTRCVE AAKDASSVKR RRLNRIPIEP HPKEEPKHPE KEQQEGVPNS DEGYWDSTTP GPEDSSSSG KKAGIPRDSY SGDALYDLIA DPDGSPATLP GGKDNEETSS LSRLKPVSPG TITCPLRTPG SLLKDSKIPI SIKHLTNLPS SHPVVHQQPS RSEMPRTKIP VSKVLVRRVS

image



Background

Alternative Names APC membrane recruitment protein 2 (Amer2) (Protein FAM123A)

Product Name: Recombinant Human AMER2 Protein
Catalog #: PEH7138



Note

For Research Use Only , Not for Diagnostic Use