


Product Name: Recombinant Human SLC26A7 Protein
Catalog #: PEH7142



Summary

Name	S26A7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TE54
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KNMKEMEFKV KTEMSETLQ QVKIISINNP LVFLNAKKFY TDLNMNIQKE NACNQPLDDI SKCEQNTLLN SLSNGNCNEE ASQSCPNEKC YLILDCSGFT FFDYSGVSMML VEVYMDCKGR SVDVLLAHCT ASLIKAMTTY GNLDSEKPIF FESVSAASH

image



Background

Alternative Names Anion exchange transporter (Solute carrier family 26 member 7)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SLC26A7 Protein
Catalog #: PEH7142

