


Product Name: Recombinant Human RAG2 Protein
Catalog #: PEH7143



Summary

Name	RAG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P55895
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKNDTIYILG GHSLANNIRP ANLYRIRVDL PLGSPAVNCT VLPGGISVSS AILTQTNNDE FVIVGGYQLE NQKRMICNII SLEDNKIEIR EMETPDWTPD IKHSKIWFSG NMGNQTVFLG IPGDNKQVVS EGFYFYMLKC AEDDTNEEQT TFTNSQTSTE DPGDSTPFED SEEFCSAEA NSFDDGDEFD TYNEDDEEDE

image



Background

Alternative Names V(D)J recombination-activating protein 2 (RAG-2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RAG2 Protein
Catalog #: PEH7143

