


Product Name: Recombinant Human TRAF3IP3 Protein
Catalog #: PEH7145



Summary

Name	T3JAM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y228
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLSTLIQASD SSWKGQLNED KLKGKLSLE NQLYTCTQKY SPWGMKKVLL EMEDQKNSYE QKAKESLQKV LEEKMNAEQQ LQSTQRLAL AEQKCEEWS QYEALKEDWR TLGTQHRELE SQLHVLQSKL QGADSRDLQM NQALRFLENE HQQQLAKIEC LQGDRDLCSL DTQDLQDQLK RSEAEKLTIV TRVQQLQGLL QNQSLQLQEQ EKLLTKKDQA LPVWSPKSF NEVEPEGTGK EKDWDLRDQL QKKTLLQLQAK EKECRELHSE LDNLSDEYLS

image



Background

Alternative Names	TRAF3-interacting JNK-activating modulator (TRAF3-interacting protein 3)
--------------------------	--

Product Name: Recombinant Human TRAF3IP3 Protein
Catalog #: PEH7145



Note

For Research Use Only , Not for Diagnostic Use