


Product Name: Recombinant Human ACYP2 Protein
Catalog #: PEH7154



Summary

Name	ACYP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P14621
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSTAQSLKSV DYEVFGRVQG VCFRMYTEDE ARKIGVVGWV KNTSKGTVTG QVQGPEDKVN SMKSWLSKVG SPSSRIDRTN FSNEKTISKL EYSNFSIRY

image



Background

Alternative Names	Acylphosphatase-2 (EC 3.6.1.7) (Acylphosphatase, muscle type isozyme) (Acylphosphate phosphohydrolase 2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use