


Product Name: Recombinant Human JMJD1C Protein
Catalog #: PEH7158



Summary

Name	JHD2C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15652
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WICEKHILWL KDYNSSNWK LFKECWKQGQ PAVVSGVHKK MNISLWKAES ISLDFGDHQA DLLNCKDSII SNANVKEFWD GFEEVSKRQK NKSGETVVLK LKDWPSGEDF KTMMPARYED LLKSLPLPEY CNPEGKFNLA SHLPGFFVRP DLGPRLCSAY GVVAAKDHDH GTTNLHIEVS DVVNILVYVG IAKGNGILSK AGILKKFEEE DLDDILRKRL KDSSEIPGAL WHIYAGKDVD KIREFLQKIS KEQGLEVLPE HDPIRDQSWY VNKKLRQRL EEYGVRTCTL IQFLGDAIVL PAGALHQVQN FHSCIQVTED FVSPEHLVES FHLTQELRL KEEINYDDKL QVKNILYHAV KEMVRALKIH EDEVEDMEEN

image



Background

Alternative Names	Probable JmjC domain-containing histone demethylation protein 2C (EC
--------------------------	--

Product Name: Recombinant Human JMJD1C Protein
Catalog #: PEH7158



1.14.11.-) (Jumonji domain-containing protein 1C) (Thyroid receptor-interacting protein 8) (TR-interacting protein 8) (TRIP-8)

Note

For Research Use Only , Not for Diagnostic Use