


**Product Name: Recombinant Human PTGFRN Protein**  
**Catalog #: PEH7163**



## Summary

|                              |                                                                                                                                                                                                                |
|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | FPRP_HUMAN                                                                                                                                                                                                     |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                            |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                               |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                          |
| <b>Accession #</b>           | Q9P2B2                                                                                                                                                                                                         |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                         |
| <b>Species</b>               | Human                                                                                                                                                                                                          |
| <b>Molecular Mass</b>        | 15 kDa                                                                                                                                                                                                         |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                      |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                             |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| <b>Sequence</b>              | VVKARQPKPF FAAGNTFEMT CKVSSKNIKS PRYSVLIMAE KPVGDLSSPN<br>ETKYIISLDQ DSVVKLENWT DASRVDGVVL EKVQEDEFY RMYQTQVSDA<br>GLYRCMVTAW                                                                                  |

## image



## Background

|                          |                                                                                                                                                                                                                                                         |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | Prostaglandin F2 receptor negative regulator (CD9 partner 1) (CD9P-1) (Glu-Trp-Ile EWI motif-containing protein F) (EWI-F) (Prostaglandin F2-alpha receptor regulatory protein) (Prostaglandin F2-alpha receptor-associated protein) (CD antigen CD315) |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## Note

**Product Name: Recombinant Human PTGFRN Protein**  
**Catalog #: PEH7163**



---

For Research Use Only , Not for Diagnostic Use