


Product Name: Recombinant Human SWAP70 Protein
Catalog #: PEH7164



Summary

Name	SWP70_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UH65
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VLKQGYMMKK GHRRKNWTER WFLKPNIIIS YYVSEDLKDK KGDILLDENC CVESLPDKDG KKCLFLVKCF DKTFEISASD KKKKQEWIQA IHSTIHLKL GSPPPHKEAR QRRKELRKKQ LAEQEELERQ MKELQAANES KQEQLEAVRK KLEEAASRAA EEEKRLQTQ VELQARFSTE LEREKLIRQQ MEEQVAQKSS ELEQYLQVRV ELEDMYLKLQ EALEDERQAR QDEETVRKLQ ARLLEEESK RAELEKWHLE QQQAIQTTEA EKQELENQRV LKEQALQEAM EQLEQLELER KQALEQYEEV KKKLEMATNK TSWKDKVAH HEGLIRLIEP GSKNPHLITN WGPAAFTEAE LEEREKNWKE

image



Background

Alternative Names	Switch-associated protein 70 (SWAP-70)
--------------------------	--

Product Name: Recombinant Human SWAP70 Protein
Catalog #: PEH7164



Note

For Research Use Only , Not for Diagnostic Use