


Product Name: Recombinant Human DAPK2 Protein
Catalog #: PEH7165



Summary

Name	DAPK2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UIK4
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MFQASMRSPN MEPFKQQKVE DFYDIGEELG SGQFAIVKKC REKSTGLEYA AKFIKKRQSR ASRRGVSREE IEREVSILRQ VLHHNVITLH DVYENRTDVV LILELVSGGE LFDFLAQKES LSEEEATSFI KQILDGVNYL HTKKIAHFDL KPENIMLLDK NIPIPHIKLI DFGLAHEIED GVEFKNIFGT PEFVAPEIVN YEPLGLEADM WSIGVITYIL LSGASPFLGD TKQETLANIT AVSYDFDEEF FSQTSELAKD FIRKLLVKET RKRLTIQEAL RHPWITPVDN QQAMVRRESV VNLENFRKQY VRRRWKLSFS IVSLCNHLTR SLMKKVHLRP DEDLRNCESD TEEDIARRKA LHPRRRSSTS

image



Background

Alternative Names	Death-associated protein kinase 2 (DAP kinase 2) (EC 2.7.11.1) (DAP-kinase-
--------------------------	---

Product Name: Recombinant Human DAPK2 Protein
Catalog #: PEH7165



related protein 1) (DRP-1)

Note

For Research Use Only , Not for Diagnostic Use