


**Product Name: Recombinant Human PPIE Protein**  
**Catalog #: PEH7167**



## Summary

<b>Name</b>	PPIE_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9UNP9
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	34 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MATTKRVLYV GGLAEEVDDK VLHAAFIPFG DITDIQIPLD YETEKHRGFA FVEFELAEDA AAAIDNMNES ELFGRTIRVN LAKPMRIKEG SSRPVWSDDD WLKKFSGKTL EENKEEEGSE PPKAETQEGE PIAKKARSNP QVYMDIKIGN KPAGRIQMLL RSDVVPMTAE NFRCLCTHEK GFGFKGSSFH RIIPQFMCQG GDFTNHNGTG GKSIYGKKFD DENFILKHTG PGLLSMANSV PNTNGSQFFL TCDKTDWLDG KHVVFGEVTE GLDVLRQIEA QGSKDGKPKQ KVIADCGEY V

## image



## Background

**Alternative Names** Peptidyl-prolyl cis-trans isomerase E (PPIase E) (EC 5.2.1.8) (Cyclophilin E) (Cyclophilin-33) (Rotamase E)

**Product Name: Recombinant Human PPIE Protein**  
**Catalog #: PEH7167**

---



**Note**

For Research Use Only , Not for Diagnostic Use