


Product Name: Recombinant Human PIK3CD Protein
Catalog #: PEH7169



Summary

Name	PK3CD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00329
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NPAPQVQKPR AKPPPIPAKK PSSVSLWSLE QPFRIELIQG SKVNADERMK LVVQAGLFHG NEMLCKTVSS SEVSVCSEPV WKQRLEFDIN ICDLPRMARL CFALYAVIEK AKKARSTKKK SKKADCPIAW ANLMLFDYKD QLKTGERCLY MWPSVPDEKG ELLNPTGTVR SNPNTDSAAA LLICLPEVAP

image



Background

Alternative Names	Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit delta isoform (PI3-kinase subunit delta) (PI3K-delta) (PI3Kdelta) (PtdIns-3-kinase subunit delta) (EC 2.7.1.153) (Phosphatidylinositol 4,5-bisphosphate 3-kinase 110 kDa catalytic subunit delta) (PtdIns-3-kinase subunit p110-delta) (p110delta)
--------------------------	---

Product Name: Recombinant Human PIK3CD Protein
Catalog #: PEH7169



Note

For Research Use Only , Not for Diagnostic Use