



Product Name: Recombinant Human CACNA1A Protein
Catalog #: PEH7170



Summary

Name	CAC1A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00555
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RNNNTNKSRA AEPTVDQRLG QQRAEDFLRK QARYHDRARD PSGSAGLDAR RPWAGSQEAE LSREGPYGRE SDHHAREGSL EQPGFWEGEA ERGKAGDPHR RHHVHRQGGSR ESRSGSPRTG ADGEHRRHRA HRRPGEEGPE DKAERRARHR EGSRPARGGE GEGEGPDGGE

image



Background

Alternative Names	Voltage-dependent P/Q-type calcium channel subunit alpha-1A (Brain calcium channel I) (BI) (Calcium channel, L type, alpha-1 polypeptide isoform 4) (Voltage-gated calcium channel subunit alpha Cav2.1)
--------------------------	--

Note

Product Name: Recombinant Human CACNA1A Protein
Catalog #: PEH7170



For Research Use Only , Not for Diagnostic Use