


Product Name: Recombinant Human FBP2 Protein
Catalog #: PEH7171



Summary

Name	F16P2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00757
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MTDRSPFETD MLTLTRYVME KGRQAKGTGE LTQLLSMLT AIKAISSAVR KAGLAHLYGI AGSVNVTGDE VKKLDVLSNS LVINMVQSSY STCVLVSEEN KDAIITAKEK RGKYVVC FDP LDGSSNIDCL ASIGTIFAIY RKTSEDEPSE KDALQCGRNI VAAGYALYGS ATLVALSTGQ GVDLFMLDPA LGEFVLVEKD VKIKKKGKIY SLNEGAYKYF DAATTEYVQK KKFPEDEGSAP YGARYVGSMV ADVHRTLVIY GIFLYPANQK SPKGKLRLLY ECNPVAYIIE QAGGLATTGT QPVLVDVKPEA IHQRVPLILG SPEDVQEYLT CVQKNQAGS

image



Background

Alternative Names	Fructose-1,6-bisphosphatase isozyme 2 (FBPase 2) (EC 3.1.3.11) (D-fructose-1,6-bisphosphate 1-phosphohydrolase 2) (Muscle FBPase)
--------------------------	---

Product Name: Recombinant Human FBP2 Protein
Catalog #: PEH7171



Note

For Research Use Only , Not for Diagnostic Use