


Product Name: Recombinant Human CACNA1G Protein
Catalog #: PEH7174



Summary

Name	CAC1G_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43497
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RGSLEREAKS SFDLPDTLQV PGLHRTASGR GSASEHQDCN GKSASGRLAR ALRPDDPPLD GDDADDEGNL SKGERVRAWI

image



Background

Alternative Names	Voltage-dependent T-type calcium channel subunit alpha-1G (Cav3.1c) (NBR13) (Voltage-gated calcium channel subunit alpha Cav3.1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use