


Product Name: Recombinant Human KCNQ2 Protein
Catalog #: PEH7176



Summary

Name	KCNQ2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43526
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DRVFSSPRGV AAKGKGSPQA QTVRRSPSAD QSLEDSPSKV PKSWSFGBDRS RARQAFRIKG AASRQNSEEA SLPGEDIVDD

image



Background

Alternative Names	Potassium voltage-gated channel subfamily KQT member 2 (KQT-like 2) (Neuroblastoma-specific potassium channel subunit alpha KvLQT2) (Voltage-gated potassium channel subunit Kv7.2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use