


Product Name: Recombinant Human HCRT2 Protein
Catalog #: PEH7177



Summary

Name	OX2R_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43614
Host	E.coli
Species	Human
Molecular Mass	7 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PCRNWSSASE LNETQEPFLN PTDYDDEEFL RYLWREYLHP

image



Background

Alternative Names Orexin receptor type 2 (Ox-2-R) (Ox2-R) (Ox2R) (Hypocretin receptor type 2)

Note

For Research Use Only , Not for Diagnostic Use