


Product Name: Recombinant Human GUCY1B2 Protein
Catalog #: PEH7180



Summary

Name	GUCYB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75343
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KVAAGEFKSC TILFSDVVF TNICTACEPI QIVNVLSMY SKFDRLTSVH AVYKVETIGD AYMVVGVPV PIGNHAQRVA NFALGMRISA KEVTNPVTGE PIQLRVGIHT GPVLADVVD KMPRYCLFGD TVNTASRMES HGLPNKVHLS PTAYRALKNQ GFKIIERGEI EVKGKGRMTT YFLIQNLNAT EDEIMGRSKT

image



Background

Alternative Names Guanylate cyclase soluble subunit beta-2 (GCS-beta-2) (EC 4.6.1.2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GUCY1B2 Protein
Catalog #: PEH7180

