


Product Name: Recombinant Human DDN Protein
Catalog #: PEH7184



Summary

Name	DEND_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O94850
Host	E.coli
Species	Human
Molecular Mass	74 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GWRESLGLGE GAGPETLEGW KATTRAHTLP RSSQGLSRGE GVFVIDATCV VIRSQYVPTP RTQQVQLLPS GVTRVVGDSP SQSKPGKEEG EGATVFPSPC QKRLSSRLL HQPGGGRGGE AEGGRP GDST LEERTFRILG LPAPEVNLRD APTQPGSPEH QALGPAASGA QGRAEGSEVA VVQRRAGR GW ARTPGPYAGA LREAVSRIRR HTAPDSDTDE AEELSVHSGS SDGSDTEAPG ASWRNERTLP EVGNSSPEED GKTAELSDSV GEILDVISQT EEVLFGVRDI RGTQQGNRKR

image



Background

Alternative Names Dendrin

Product Name: Recombinant Human DDN Protein
Catalog #: PEH7184



Note

For Research Use Only , Not for Diagnostic Use