


Product Name: Recombinant Human TOX Protein
Catalog #: PEH7185



Summary

Name	TOX_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94900
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QPAAAPDAPC LGPSPCLDPY YCNKFDGENM YMSMTEPSQD YVPASQSYPG PSLESEDFNI PPITPPSLPD HSLVHLNEVE SGYHSLCHPM NHNGLLPFHP QNMDLPEITV SNMLGQDGTI LSNSISVMPD IRNPEGTQYS SHPQMAAMRP RGQPADIRQQ PGMMPHGQLT TINQSQLSAQ LGLNMGGSNV PHNSPSPPGS KSATPSPSSS VHEDEGDDTS KINGGEKRPA SDMGKKPKTP KKKKKKDPNE PQKPV SAYAL FFRDTQAAIK GQNP NATFGE VSKIVASMWD GLGEEQKQVY KKKTEAAKKE

image



Background

Alternative Names Thymocyte selection-associated high mobility group box protein TOX (Thymus high mobility group box protein TOX)

Product Name: Recombinant Human TOX Protein
Catalog #: PEH7185



Note

For Research Use Only , Not for Diagnostic Use