


**Product Name: Recombinant Human NFAT5 Protein**  
**Catalog #: PEH7186**



## Summary

<b>Name</b>	NFAT5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	O94916
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	32 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	GQYPVKSE GK ELKIVVQPET QHRARYLTEG SRGSVKDRTQ QGFPTVKLEG HNEPVVLQVF VGNDSGRVKP HGFYQACRVT GRNTTPCKEV DIEGTTVIEV GLDPSNNMTL AVDCVGILKL RNADVEARIG IAGSKKKSTR ARLVFRVNIM RKDGSTLTLQ TPSSPILCTQ PAGVPEILKK SLHSCSVKGE EEVFLIGKNF LKGTKVIFQE NVSDENSWKS EAEIDMELFH QNHLIVKVPP YHDQHITLPV SVGIVVTNA GRSHDVQPFT

## image



## Background

<b>Alternative Names</b>	Nuclear factor of activated T-cells 5 (NF-AT5) (T-cell transcription factor NFAT5) (Tonicity-responsive enhancer-binding protein) (TonE-binding protein) (TonEBP)
--------------------------	---

**Product Name: Recombinant Human NFAT5 Protein**  
**Catalog #: PEH7186**

---



**Note**

For Research Use Only , Not for Diagnostic Use