


Product Name: Recombinant Human IPO7 Protein
Catalog #: PEH7189



Summary

Name	IPO7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O95373
Host	E.coli
Species	Human
Molecular Mass	104 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DQNLQTALEL TRRCLIDRE MPVKVEAAIA LQVLISNQEK AKEYITPFIR PVMQALLHII RETENDDLTN VIQKMICEYS EEVTPIAVEM TQHLAMTFNQ VIQTGPDEEG SDDKAVTAMG ILNTIDTLLS VVEDHKEITQ QLEGICLQVI GTVLQQHVLE FYEEIFSLAH SLTCQQVSPQ MWQLLPLVFE VFQQDGFDFY TDMMPLLHNY VTVDTDTLLS DTKYLEMIYS MCKKVLTVGA GEDAECHAAK LLEVILQCK GRGIDQCIPL FVEAALERLT REVKTSELRT MCLQVAIAAL YYNPHLLNT LENLRFNNV EPVTNHFITQ WLNDVDCFLG LHDRKMCVLG LCALIDMEQI PQVLNQVSGQ ILPAFILLFN GLKRAYACHA EHENDSDDDD EAEDDDETEE LGSEDDIDE DGQEYLEILA KQAGEDGDDE DWEEDDAEET ALEGYSTIID DEDNPVDEYQ IFKAIFQTIQ NRPVWYQAL THGLNEEQRK QLQDIATLAD QRRAAHESKM IEKHGGYKFS APVVPSSFNF

image



Background

Product Name: Recombinant Human IPO7 Protein
Catalog #: PEH7189



Alternative Names

Importin-7 (Imp7) (Ran-binding protein 7) (RanBP7)

Note

For Research Use Only , Not for Diagnostic Use