


Product Name: Recombinant Human GLP2R Protein
Catalog #: PEH7192



Summary

Name	GLP2R_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95838
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FRFLGKCPKK LSEGDAEKL RKLQPSLNSG RLLHLAMRGL GELGAQPQQD HARWPRGSSL SECSEGDVTM ANTMEEILEE

image



Background

Alternative Names Glucagon-like peptide 2 receptor (GLP-2 receptor) (GLP-2-R) (GLP-2R)

Note

For Research Use Only , Not for Diagnostic Use