


Summary

Name	GSHR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P00390
Host	E.coli
Species	Human
Molecular Mass	53 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MACRQEPQPQ GPPPAAGAVA SYDYLVIIGGG SGGLASARRA AELGARAADV ESHKLGSTCV NVGCVPKKVM WNTAVHSEFM HDHADYGFPS CEGKFNWRVI KEKRDAYVSR LNAIQNNLT KSHIEIRGH AAFTSDPKPT IEVSGKKYTA PHILIATGGM PSTPHESQIP GASLGITSDG FFQLEELPGR SVIVGAGYIA VEMAGILSAL GSKTSLMIRH DKVLRSFDSM ISTNCTEELE NAGVEVLKFS QVKEVKKTLS GLEVSMVTAV PGRLPVMTMI PDVDCLLWAI GRVPNTKDLS LNKLGITDD KGHIIVDEFQ NNTNVKGIYAV GDVCGKALLT PVAIAAGRKL AHRLEFYKED SKLDYNNIPT VVFSHPPIGT VGLTEDEAIH KYGIENVKTY STSFPMYHA VTKRKTCKVM KMVCANKEEK VVGIHMQGLG CDEMLQGFV AVKMGATKAD FDNTVAIHPT SSEELVTLR

image



Background

Product Name: Recombinant Human GSR Protein
Catalog #: PEH7194



Alternative Names

Glutathione reductase, mitochondrial (GR) (GRase) (EC 1.8.1.7)

Note

For Research Use Only , Not for Diagnostic Use