


Product Name: Recombinant Human CHR2 Protein
Catalog #: PEH7210



Summary

Name	ACM2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08172
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RASKSRIKKD KKEPVANQDP VSPSLVQGRI VKPNNNNMPS SDDGLEHNKI QNGKAPRPV TENCVQGEEK ESSNDSTSVS AVASNMRDDE ITQDENTVST SLGHSKDENS KQTCIRIGTK TPKSDSCTPT NTTVEVVGSS GQNGDEKQNI VARKIVKMTK QPAKKKPPPS

image



Background

Alternative Names Muscarinic acetylcholine receptor M2

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CHRM2 Protein
Catalog #: PEH7210

