



Product Name: Recombinant Human ADRB1 Protein
Catalog #: PEH7211



Summary

Name	ADRB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08588
Host	E.coli
Species	Human
Molecular Mass	8 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AQKQVKKIDS CERRFLGGPA RPPSPSPSPV PAPAPPPGPP RPAAAAATAP LANGRAGKRR

image



Background

Alternative Names Beta-1 adrenergic receptor (Beta-1 adrenoreceptor) (Beta-1 adrenoceptor)

Note

For Research Use Only , Not for Diagnostic Use