


Product Name: Recombinant Human ADRA2A Protein
Catalog #: PEH7212



Summary

Name	ADA2A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P08913
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FAPCLIMILV YVRIYQIAKR RTRVPPSRRG PDAVAAPPGG TERRPNGLGP ERSAGPGGAE AEPLPTQLNG APGEPAPAGP RDTDALDLEE SSSSDHAERP PGRRRPERGP RGKGKARASQ VKPGDSLPRR GPGATGIGTP AAGPGEERVG

image



Background

Alternative Names	Alpha-2A adrenergic receptor (Alpha-2 adrenergic receptor subtype C10) (Alpha-2A adrenoreceptor) (Alpha-2A adrenoceptor) (Alpha-2AAR)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use