


Product Name: Recombinant Human ODC1 Protein
Catalog #: PEH7218



Summary

Name	DCOR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P11926
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MNNFGNEEFD CHFLDEGFTA KDILDQKINE VSSDDKDAF YVADLGDLK KHLRWLKALP RVTPFYAVKC NDKAIVKTL AATGTGFDCA SKTEIQLVQS LGVPPERIIY ANPCKQVSQI KYAANNGVQM MTFDSEVELM KVARAHPKAK LVLRIATDDS KAVCRLSVKF GATLRTRSLL LERAKELNID VVGVSFHVGS GCTDPETFVQ AISDARCVFD MGAEVGFMSY LLDIGGGFPG SEDVKLKFEF ITGVINPALD KYFPSDSGVR IIAEPGRYYV ASFTLAVNI IAKKIVLKEQ TGSDDEDESS EQTFMYYVND GYGSFNCIL YDHAHVKPLL QKRPKPDEKY YSSSIWGPTC DGLDRIVERC DLPEMHVGDW MLFENMGAYT VAAASTFNFG QRPTIYYVMS GPAWQLMQQF QNPDPPEVE EQDASTLPVS CAWESGMKRH RAACASASIN V

image



Background

Product Name: Recombinant Human ODC1 Protein
Catalog #: PEH7218



Alternative Names

Ornithine decarboxylase (ODC) (EC 4.1.1.17)

Note

For Research Use Only , Not for Diagnostic Use