


**Product Name: Recombinant Human MYH8 Protein**  
**Catalog #: PEH7221**



## Summary

<b>Name</b>	MYH8_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P13535
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	29 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	SRKIAEQELL DASERVQLLH TQNTSLINTK KKLENDVSQL QSEVEEVIQE SRNAEEKAKK AITDAAMMAE ELKKEQD TSA HLERMKNLE QTVKDLQHRL DEAEQLALKG GKKQIQKLEA RVRELEGEVE NEQKRNAEAV KGLRKHERRV KELTYQTEED RKNVLRQLDL VDKLQAKVKS YKRQAEAEAE QSNANLSKFR KLQHELEEAE ERADIAESQV NKLRVKSREV

## image



## Background

**Alternative Names** Myosin-8 (Myosin heavy chain 8) (Myosin heavy chain, skeletal muscle, perinatal) (MyHC-perinatal)

## Note

**Product Name: Recombinant Human MYH8 Protein**  
**Catalog #: PEH7221**



---

For Research Use Only , Not for Diagnostic Use