



Product Name: Recombinant Human CHR3 Protein
Catalog #: PEH7227



Summary

Name	ACM3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P20309
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RTKELAGLQA SGTEAETENF VHPTGSSRSC SSYELQQSM KRSNRRKYGR CHFVFTTKSW KPSSEQMDQD HSSDSWNNN DAAASLENSA SSDEEDIGSE TRAIYSIVLK LPGHSTILNS TKLPSSDNLQ VPEELGMVD LERKADKLQA QKSVDDGGSF PKSFSKLIQ LESAVDTAKT SDVNSSVGKS TATLPLSFKE ATLAKRFALK TRSQITKRKR MSLVKEKAA

image



Background

Alternative Names Muscarinic acetylcholine receptor M3

Note

Product Name: Recombinant Human CHRM3 Protein
Catalog #: PEH7227



For Research Use Only , Not for Diagnostic Use