


Product Name: Recombinant Human FLNA Protein
Catalog #: PEH7228



Summary

Name	FLNA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P21333
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SAAGAAPGGG VDTRDAEMPA TEKDLAEDAP WKKIQNTFT RWCNEHLKCV SKRIANLQTD LSDGLRLIAL LEVLSQKKMH RKHNQRPTFR QMQLENVSPA LEFLDRESIK LVSIDSKAIV DGNLKLILGL IWTILHYSI SMPMWDEEED EEAKKQTPKQ RLLGWIQNKL PQLPITNFSR DWQSGRALGA LVDSCAPGLC PDWDSWDASK PVTNAREAMQ QADDWLGIPIQ VITPEEIVDP NVDEHSVMTY LSQFPKAKLK

image



Background

Alternative Names	Filamin-A (FLN-A) (Actin-binding protein 280) (ABP-280) (Alpha-filamin) (Endothelial actin-binding protein) (Filamin-1) (Non-muscle filamin)
--------------------------	--

Product Name: Recombinant Human FLNA Protein
Catalog #: PEH7228



Note

For Research Use Only , Not for Diagnostic Use