


Product Name: Recombinant Human FLNA Protein
Catalog #: PEH7229



Summary

Name	FLNA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P21333
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FQVDTSKAGV APLQVKVQGP KGLVEPVDVV DNADGTQTVN YVPSREGPYS ISVLYGDEEV PRSPFKVKVL PTHDASKVKA SGPLNNTTGV PASLPVEFTI DAKDAGEGLL AVQITDPEGK PKKTHIQDNH

image



Background

Alternative Names Filamin-A (FLN-A) (Actin-binding protein 280) (ABP-280) (Alpha-filamin) (Endothelial actin-binding protein) (Filamin-1) (Non-muscle filamin)

Note

For Research Use Only , Not for Diagnostic Use