


Product Name: Recombinant Human ACAT1 Protein
Catalog #: PEH7233



Summary

Name	THIL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P24752
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VSKPTLKEVV IVSATRTPIG SFLGSLLLP ATKLGSIQ GAIEKAGIPK EEVKEAYMGN VLQGGEGQAP TRQAVLGAGL PISTPCTTIN KVCASGMKAI MMASQSLMCG HQDVMVAGGM ESMSNVPYVM NRGSTPYGGV KLEDLIVKDG LTDVYNKIHM GSCAENTAKK LNIARNEQDA YAINSYTRSK AAWEAGKFGN EVIPVTVTVK GQPDVVVKED EEYKRVDFSK VPKLKTVFQK ENGTVTAANA STLNDGAAAL VLMTADAAKR LNVTPLARIV AFADAAVEPI DFPIAPVYAA SMVLKDVGLK KEDIAMWEVN EAFSLVLAN IKMLEIDPQK VNINGGAVSL GHPIGMSGAR IVGHLTHALK QGEYGLASIC NGGGGASAML IQKL

image



Background

Alternative Names Acetyl-CoA acetyltransferase, mitochondrial (EC 2.3.1.9) (Acetoacetyl-CoA)

Product Name: Recombinant Human ACAT1 Protein
Catalog #: PEH7233



thiolase) (T2)

Note

For Research Use Only , Not for Diagnostic Use