


Product Name: Recombinant Human WEE1 Protein
Catalog #: PEH7238



Summary

Name	WEE1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P30291
Host	E.coli
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDLLPGACP GADEAGGGAE GDSWEEEGFG SSSPVKSPAA PYFLGSSFSP VRCGGPGDAS PRGCGARRAG EGRRSPRPDH PGTPPHKTFR KLRLFDTPHT PKSLLSKARG IDSSSVKLRG SSLFMDTEKS GKREFDVRQT PQVNINPFTP DSLLHSSGQ CRRRKRTYWN DSCGEDMEAS DYELEDETRP AKRITITESN MKSRYTTEFH ELEKIGSGEF GSVFKCVKRL DGCIYAIKRS KKPLAGSVDE QNALREVVYAH AVLGQHSHVV RYFSAWAEDD HMLIQNEYCN GGSLADAISE NYRIMSYFKE AELKDLLLQV GRGLRYIHSM SLVHMDIKPS NIFISRTSIP NAASEEGDED DWASNKVMFK IGD LGHVTRI SSPQVEEGDS RFLANEVLQE NYTHLPKADI FALALTVCVA AGAEPLPRNG DQWHEIRQGR LPRIPQVLSQ EFTELLKVM IHPDPERRPSA

image



Background

Product Name: Recombinant Human WEE1 Protein
Catalog #: PEH7238



Alternative Names

Wee1-like protein kinase (WEE1hu) (EC 2.7.10.2) (Wee1A kinase)

Note

For Research Use Only , Not for Diagnostic Use