


Product Name: Recombinant Human HPD Protein
Catalog #: PEH7241



Summary

Name	HPPD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P32754
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MTTYSDKGAK PERGRFLHFH SVTFWVGNAK QAASFYCSKM GFEPLAYRGL ETGSREVVSH VIKQGKIVFV LSSALNPWNK EMGDHLVKHG DGVKDIAFEV EDCDYIVQKA RERGAKIMRE PWVEQDKFGK VKFAVLQTYG DTTHTLVEKM NYIGQFLPGY EAPAFMDPLL PKLPKCSLEM IDHIVGNQPD QEMVSASEWY LKNLQFHRFW SVDDTQVHTE YSSLRSIVVA NYEESIKMPI NEPAPGKKKS QIQEYVDYNG GAGVQHIALK TEDIITAIRH LRERGLEFLS VPSTYYKQLR EKLKTAKIKV KENIDALEEL KILVDYDEKG YLLQIFTKPV QDRPTLFLEV IQRHNNHQGFG AGNFNSLFKA FEEEQNLRGN LTNMETNGVV PGM

image



Background

Alternative Names 4-hydroxyphenylpyruvate dioxygenase (EC 1.13.11.27) (4-

Product Name: Recombinant Human HPD Protein
Catalog #: PEH7241



hydroxyphenylpyruvic acid oxidase) (4HPPD) (HPD) (HPPDase)

Note

For Research Use Only , Not for Diagnostic Use