


Product Name: Recombinant Human SCN4A Protein
Catalog #: PEH7245



Summary

Name	SCN4A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35499
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FNVATEESSE PLGEDDFEMF YETWEKFDPD ATQFIAYSRL SDFVDTLQEP LRIAKPNKIK LITLDLPMVP GDKIHCLDIL FALTKEVLGD SGEMDALKQT MEEKFMAANP SKVSYEPITT TLKRKHEEVC AIKIQRAYRR HLLQRSMKQA SYMVRHSHDG SGDDAPEKEG LLANTMSKMY GHENGNSSSP SPEEKGEAGD AGPTMGLMPI SPSDTAWPPA PPPGQTVRPG

image



Background

Alternative Names	Sodium channel protein type 4 subunit alpha (SkM1) (Sodium channel protein skeletal muscle subunit alpha) (Sodium channel protein type IV subunit alpha) (Voltage-gated sodium channel subunit alpha Nav1.4)
--------------------------	--

Product Name: Recombinant Human SCN4A Protein
Catalog #: PEH7245



Note

For Research Use Only , Not for Diagnostic Use