



Product Name: Recombinant Human MYH11 Protein
Catalog #: PEH7250



Summary

Name	MYH11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35749
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLEGDLKDLE LQADSAIKGR EEAIKQLRKL QAQMKDFQRE LEDARASRDE IFATAKENEK KAKSLEADLM QLQEDLAAAE RARKQADLEK EELAEELASS LSGRNALQDE KRRLEARIAQ LEEEEEEQGG NMEAMSDRVR KATQQAEQLS NELATERSTA QKNESARQQL ERQNKELRSK LHEMEGAVKS KFKSTIAALE AKIAQLEEQV EQEAREKQAA TKSLKQKDKK LKEILLQVED ERKMAEQYKE QAEKGNARVK QLKRRQLEEAE EESQRINANR RKLQRELDEA TESNEAMGRE VNALKSKLRR GNETSFVPSR RSGGRRVIEN ADGSEETDT RDAFNGTKA

image



Background

Alternative Names	Myosin-11 (Myosin heavy chain 11) (Myosin heavy chain, smooth muscle isoform) (SMMHC)
--------------------------	---

Product Name: Recombinant Human MYH11 Protein
Catalog #: PEH7250



Note

For Research Use Only , Not for Diagnostic Use