


Product Name: Recombinant Human CTSK Protein
Catalog #: PEH7258



Summary

Name	CATK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P43235
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QLKKKTGKLL NLSPQNLVDC VSENDGCGGG YMTNAFYVQ KNRGIDSEDA YPYVGQEESC MYNPTGKAAK CRGYREIPEG NEKALKRAVA

image



Background

Alternative Names Cathepsin K (EC 3.4.22.38) (Cathepsin O) (Cathepsin O2) (Cathepsin X)

Note

For Research Use Only , Not for Diagnostic Use